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
1. Prepare JC-1 working solution:

1.1 Prepare a 2 to 10 mM stock solution of JC-1 in high-quality, anhydrous DMSO. The stock solution should be used promptly; any remaining solution should be aliquoted and frozen at ≤ -20 °C.

Note: Avoid repeated freeze-thaw cycles, and protect from light.

1.2 Prepare a 1X JC-1 working solution: On the day of the experiment, either dissolve JC-1 solid in DMSO or thaw an aliquot of the JC-1 stock solution to room temperature. Prepare a 10 to 30 μ M 1X JC-1 working solution in Hanks and 20 mM Hepes buffer (HHBS) or buffer of your choice, pH 7 with 0.02% Pluronic® F-127(MA0061). Mix them well by vortexing.

Note: JC-1 is not water soluble, so it intends to aggregate in solution. It is recommended to filter the JC-1 working solution before loading it into the cells.

2. Run JC-1 assay with a fluorescence microplate reader:

2.1 Treat cells with test compounds for a desired period of time (For example, Jurkat cells can be treated with camptothecin for 4-6 hours) to induce apoptosis. For blank wells (medium without the cells), add the corresponding amount of compound buffer.

2.2 Add 100 μ L/well/96-well plate or 25 μ L/well/384-well plate of JC-1 working solution (from Step 1.2) into the cell plate.

2.3 Incubate the JC-1 loading plate in a 37 °C, 5% CO₂ incubator for 15-60 min.

Note: The appropriate incubation time depends on the individual cell type and cell concentration used.

Optimize the incubation time for each experiment.

2.4 Remove the JC-1 working solution from the plate, wash the cells with HHBS or buffer of your choice. Add 100 μ L/well/96-well plate or 25 μ L/well/384-well plate of HHBS back to the cell plate.

2.5 Monitor the fluorescence change at Ex/Em = 490/525 nm and 490/590 nm for ratio analysis.

3. Run JC-1 assay with a fluorescence microscope or a flow cytometer:

3.1 Treat cells with test compounds for a desired period of time (For example, Jurkat cells can be treated with camptothecin for 4-6 hours) to induce apoptosis.

3.2 Centrifuge the cells to get 1-5 $\times 10^5$ cells per tube.

3.3 Resuspend cells in 500 μ L of JC-1 working solution (from Step 1.2).

3.4 Incubate at room temperature or 37° C for 10 to 30 min, protected from light.

3.5 Wash the cells with HHBS or buffer of your choice. Resuspend cells in 500 μ L of HHBS to get 1-5 $\times 10^5$ cells per tube.

3.6 Monitor the fluorescence change at Ex/Em = 490/525 nm and 490/590 nm with a fluorescence microscope (using FITC and TRITC filters) or a flow cytometer (using FL1 and FL2 channels).

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